

## Cleaning Code : **WS / Bleach Cleanable**

Water based cleaning agents and foam may be used for cleaning. This fabric may also be cleaned with mild, water-free solvents.

Cleaning by a professional service is recommended.

May be cleaned with bleach following the proper ratio of 90% water : 10% bleach.

All cleaning should allow for a 10-minute maximum sit time and then **MUST BE PROMPTLY AND PROPERLY RINSED.**

or

## Cleaning Instructions: **CRYPTON®**

Most new stains come off with a damp cloth and water. There is one suggested cleaning method that is universally and easily available - the use of enzyme detergent (listed below).

### Spot-Cleaning:

- Mix five parts water with one part enzyme laundry detergent (e.g., Tide® or Cheer®) to form a mixture that works well on most stains. Follow these simple steps:
- Remove excess soil and apply the soap mixture to the stained area
- Agitate stain with a soft brush
- Allow solution to remain on stain for one minute
- Blot up the stain with a dry, white towel, rinse any remaining soap and blot the area again
- Repeat if necessary

### If Extraction is Required:

- Before wet-cleaning, vacuum the upholstery
- Apply cleaning solution to specific spots and lightly brush
- Follow directions on filling extraction machine with proper solution, wet the area with pre-spray and vacuum up the solution
- Be sure to extract the entire surface of the upholstered furniture to achieve uniform cleanliness

### If Beach-Cleaning is Required:

- Using bleach on Crypton (as any other textile) is dependent on the fiber content. We recommend cleaning with a 10% bleach/water mixture only on solution-dyed polyester, acrylic or 100% polyester
- Spray the 10% bleach/water mixture on the fabric and allow to linger for the time specified by the Centers for Disease Control and Prevention. Make certain to thoroughly rinse after bleach-cleaning



Arc-Com Fabrics, Inc.  
33 Ramland South  
Orangeburg, NY 10962

Tel (845) 365-1100 (800) 223-5466  
Fax (845) 365-1285 [www.arc-com.com](http://www.arc-com.com)

## Cleaning Instructions: **CRYPTON<sup>®</sup>**

We know from experience that every facility has preferred cleaners and cleaning practices. Because of this fact, we've tested hundreds of products to ensure they won't harm textiles, including our Crypton fabric. If there is a special cleaning agent you plan to use with a Crypton product, or for more information, contact us at [inquiries@arc-com.com](mailto:inquiries@arc-com.com)

The following products may be used to safely clean Crypton fabric. Follow the manufacturers recommended instructions. As with all cleaners, the fabric must be rinsed thoroughly.

<b>Product Name</b>	<b>Producer</b>
#700 Disinfectant Cleaner, Quat 128	Waxie
3M Brand HB Quat Disinfectant Cleaner-RTU	3M
3M Neutral Quat Disinfectant Cleaner	3M
3M Quat Disinfectant Cleaner Concentrate	3M
3M TM HB Quat Disinfectant Concentrate	3M
A-456-N	Eco-Lab
Buckeye Sanicare Quat 256	Buckeye
Carpet & Textile Stain Digester	Sysco Chemical
Eco-Star Detergent MP	Eco-Lab
ED - Everyday Disinfectant	Wood Wyant, Inc.
Enverros Fiberfresh Ultra Pro	Eco-Lab
QC 31	Eco-Lab
Quik Fill 120	Professional Products
Super HDQ Neutral	Spartan Chemical
Sureblend Hospital Disinfectant	Amsan
Ultraquat Germidical Detergent	Wood Wyant, Inc.
Unbelievable! Stain Remover & Deodorizer	Core Products
Virex II 256 Concentrate	Johnson Diversey
Virex II 256 RTU One Step Disinfectant	Johnson Diversey
Viva	Eco-Lab